• For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website www.dreve.dentamid.de.

• Post curing

PCU LED N₂
Rapid Shape D,DII / Asiga MAX UV 6 minutes, 80 %
light power (vacuum)

FotoDent® flash
Rapid Shape D, DII 8 minutes + 1 minute gas pre-flow
Asiga MAX UV 10 minutes + 1 minute gas pre-flow

RS Cure
According to respective parameter set for FotoDent® cast.

• Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® cast. It is recommended to regularly take out the container of the production unit, to homogenise and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® cast.

When embedding FotoDent® cast, please note the manufacturer's data of the used embedding materials. To ensure that the resin combusts without any residue, a holding time of 60 minutes at 800–900 °C is recommended.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® cast material.

Possible dangers

Please note the safety instructions for processing FotoDent® cast.

Hazard statements: H317 May cause an allergic skin reaction.

Precautionary statements: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves. **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations. **P312** Call a POISON CENTER/doctor if you feel unwell.

Contains: Diurethandimethacrylate, triethylene glycol dimethacrylate, phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® cast packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Intended use

FotoDent® cast is a material based on (meth-)acrylate resin for DLP systems with 385 nm / 405 nm LED for the manufacturing of profiles that burn out without any residue for the dental casting technology. For post curing of FotoDent® cast built parts, the light curing units PCU LED N₂, FotoDent® flash und RS Cure are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® cast we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® cast should be well shaken approx. 1 hour before use.
- Please see further steps within printer manual.
- Carefully pour FotoDent® cast into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:
Rapid Shape D and DII printer 23 ± 3 °C
Asiga MAX UV 35 ± 3 °C
- Select the product specific parameter settings for FotoDent® cast for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2 x 6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2 x 6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.

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