Dreve Design Rev.070319-3934

• For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website [www.dreve.dentamid.de](http://www.dreve.dentamid.de).

• Post-curing:

PCU LED N <sub>2</sub> Rapid Shape D, DII	10 minutes, 80 % light output (incl. nitrogen)
Asiga MAX UV	8 minutes, 80 % light output (incl. nitrogen)

FotoDent® flash Rapid Shape D, DII	12 minutes + 1 minute gas pre-flow
Asiga MAX UV	10 minutes + 1 minute gas pre-flow

• Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® tray. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® tray.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® tray material.

### Contra indication

Intended use of FotoDent guide is manufacturing dental drill guides, an allergic reaction might occur in case of allergy to (meth)acrylates.

### Possible dangers

Please note the safety instructions for processing FotoDent® tray.

**Hazard statements: H318** Causes serious eye damage. **H317** May cause an allergic skin reaction. **H412** Harmful to aquatic life with long lasting effects.

**Precautionary statements: P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P310** Immediately call a POISON CENTER/doctor. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

**Contains:** 2-Propenoic acid, reaction product with Pentaerythrite, 1,4 Butanediol dimethacrylate, phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide, Pigments

### Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® tray packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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## Intended use

FotoDent® tray is a material based on (meth-)acrylate resin for DLP systems with 385 nm / 405 nm LED for the production of dental impression trays. For the post curing of FotoDent® tray built parts, the light curing units PCU LED N<sub>2</sub> and FotoDent® flash are recommended, which are optimal balanced to the production process.

## Indications for use

For manufacturing a custom impression tray for applying impression material.

## Processing

- During the application of FotoDent® tray we recommend wearing personal protective equipment e.g. suitable gloves, safety goggles etc.
- FotoDent® tray should be well shaken approx. 1 hour before use.
- Please see further steps within printer manual.
- Carefully pour FotoDent® tray into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:
 

Rapid Shape D and DII printer	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C
- Select the product specific parameter settings for FotoDent® tray for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2x6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.

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**Dreve Dentamid GmbH**