I. Introduction

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
I. Introduction

OPERA SYSTEM Viva HT is a pre-sintered zirconia disc 98.5 mm in diameter, which contains a plastic ring. This is designed for all milling systems using this generic-type disc. (Please refer to your milling system's technical instructions for correct machine operation.) OPERA SYSTEM Viva HT has five available thicknesses: (10mm, 14mm, 18mm, 22mm and 26mm) and has three shade variations: (HT0W, HT2W and HT3W). OPERA SYSTEM Viva HT is higher in translucency. OPERA SYSTEM Viva HT is recommended for use in fabricating FCZ (Full Contour Zirconia) restorations or the frameworks.


II. Intended Use

OPERA SYSTEM Viva HT is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, FCZ bridges, inlays, onlays and veneers).

III. Sintering Program


Sintering Program 1

Temperature	Programming Rate	Holding Time
Room Temp. — 1500°C (2732°F)	10°C/min. (18°F/min.)	—
1500°C (2732°F)	—	2 hrs
1500°C (2732°F) — Room Temp.	-10°C/min. (-18°F/min.)	—

Sintering Program 2

Temperature	Programming Rate	Holding Time
Room Temp. — 1515°C (2759°F)	35°C/min. (63°F/min.)	—
1515°C (2759°F)	—	0.5 hrs
1515°C (2759°F) — Room Temp. (*)	-45°C/min. (-81°F/min.)	—

IV. Composition

ZrO₂, Y₂O₃ etc.

V. Type and Class (ISO6872:2015)

Type:II/ Class:5

VI. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 9,9x10⁻⁶/K

VII. Directions for Use

- Take the disc from the packaging and confirm that the disc does not have a crack or other damage.
- Place the disc into the milling machine; then begin the milling process following the milling systems technical instructions.
- After milling, remove the restorations from the disc with a diamond bur, etc.
- The cutting waste or dust, which is attached to the restorations, can be removed with a gentle air stream.
- Put the restorations into the refractory sagger tray and place them in the sintering furnace.
- Depending on the performance of the sintering furnace used, review the sintering schedule shown above (III. Sintering Program) before sintering the restorations.
- After sintering, adjust the restorations with a diamond bur as needed.
- Confirm that the restorations have no cracks.
- 9)-1 FCZ:
 - Glaze baking: Create a high shine surface by polishing, especially on the contact areas, then apply the glaze on all surfaces in the usual manner.
- 9)-2 Frameworks: Build-up the dental porcelain (CERABIEN ZR or CZR PRESS LF) on the frameworks following the manufacturer's technical instructions. Check the coefficient of thermal expansion of the porcelain in the manufacturer's technical instructions to confirm compatibility.

VIII. Remarks on Handling

Contraindications:

- If the patient is hypersensitive to zirconia or any other components, this product must not be used.

Warning:

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.

Caution:

- This product should not be used when malocclusion, clenching or bruxism conditions are present.
- When milling the disc or cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the disc or cutting, grinding and polishing the restorations, use safety glasses to prevent dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
- Do not use for any purposes except for dental restoration. This product is for dental application only.
- Do not touch the items heated in the furnace with your bare hands.
- Do not take the restorations out of the sintering furnace during high temperature, as the quenching causes the breaking. However, if an elevating type furnace is used under sintering program 2, the restorations may be removed from the furnace at 800°C (1472°F) or less. When removed from the furnace, the restorations must be put on a tray made of ceramic fiber (for example: Noritake Porcelain Mat) to cool slowly.

7. There may be a large difference between the setting temperature in the program and the temperature in the real furnace when using the furnace for sintering program 2. Please use the furnace after checking with the manufacturer that the furnace and temperature listed in sintering program 2 is compatible.

8. **Sintering program 2 is only recommended for frameworks (up to 3 units), FCZ crowns, FCZ bridges (up to 3 units), inlays, onlays and veneers.**
9. The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.
10. When preparing teeth, avoid the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercuts, guide grooves, the formation of retentive holes, and sharp corners.
11. Keep the following thickness of this product for fabricating prosthetics:

Location & indication	Wall thickness
Anterior crown or bridge	0.4 mm or more
Veneer	0.4 mm or more
Posterior crown or bridge	0.5 mm or more
Inlay or onlay	0.5 mm or more

12. Use the following cross-sectional areas for connectors when fabricating bridges.

Location & indication	Connector cross section	
Anterior	2- or 3-unit bridges	7 mm ² or more
	more than 3-unit bridges	9 mm ² or more
Posterior	2- or 3-unit bridges	9 mm ² or more
	more than 3-unit bridges	9 mm ² or more

13. Keep the connected bridge pontics maximum 2 teeth in bridge constructions. When the pontic of 2 teeth continues, keep the connector cross section between the pontic at 12 mm².
14. Keep the cantilever bridge to 1 pontic tooth and the connector cross section at 12 mm².
15. Choose a shade color that is brighter than the intended color for a thick restoration, as it may look duller depending on the thickness of the restorations.
16. When using a sintering furnace for the first time and changing a sintering condition, colors after sintering may vary. Sinter a small piece of Zirconia beforehand and confirm the color.

Caution in conjunction with usage:

- Do not use the disc if there is a crack noticed after removing it from the package.
- Dry milling is recommended. If wet milling is used, the translucency of the disc may be reduced.
- If you find a crack in the restorations, do not use.
- The above sintering recommendation is only a guideline; some adjustments may be required depending on each individual furnace.
- When milling the disc, use caution when approaching the milling bur with the plastic ring to prevent detachment of the disc. Mill the disc as if leaving a zirconia part contacting the internal side of the plastic ring at 2 mm.
- Cut and remove the plastic ring prior to sintering the whole disc.
- Dispose of this product as a medical waste to prevent infection.

Storage:

- Store in a cool and dry place. Keep away from direct sunlight.
- The product should be stored at 10-30°C (50-86°F).
- Do not remove the disc from its packaging during storage.
- The disc is fragile, and requires care when handling.
- The product must be stored in an appropriate place where only dental personnel have access.
- The product must be used by the expiration date indicated on the package.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

"OPERA SYSTEM" is a trademark of EUROMAX MONACO.

"CERABIEN" is a trademark of NORITAKE CO., LIMITED.

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Kuraray Europe GmbH