## INSTRUCTIONS FOR USE

# IBT Flex Resin

<b>EN</b>	Instructions for Use	<b>GA</b>	Treoracha Úsáide
<b>DE</b>	Gebrauchsanweisung	<b>RO</b>	Instrucțiuni de utilizare
<b>FR</b>	Mode d'emploi	<b>HR</b>	Upute za upotrebu
<b>IT</b>	Istruzioni per l'uso	<b>MT</b>	Istruzzjonijiet għall-Użu
<b>ES</b>	Instrucciones de uso	<b>LV</b>	Lietošanas instrukcija
<b>PT</b>	Instruções de uso	<b>LT</b>	Naudojimo instrukcijos
<b>NL</b>	Gebruiksaanwijzing	<b>ET</b>	Kasutamise juhised
<b>DA</b>	Brugsanvisning	<b>SK</b>	Návod na použitie
<b>SV</b>	Bruksanvisning	<b>SL</b>	Navodila za uporabo
<b>FI</b>	Käyttöohjeet	<b>PL</b>	Instrukcja użytkowania
<b>CS</b>	Návod k použití	<b>BG</b>	Инструкции за употреба
<b>HU</b>	Használati útmutató	<b>JP</b>	使用説明書
<b>GR</b>	Οδηγίες χρήσης	<b>ZH</b>	使用说明

Prepared: 07/24/2023 PRNT-0109 Rev 00

# INSTRUCTIONS FOR USE

# IBT Flex Resin

Keep away from sunlight

Use-by date

Consult Instructions for Use

European Conformity

Batch Code

Catalog Number

Manufacturer

Temperature Limit

Authorized representative  
in the European Community

Medical Device

Environmental hazard

Warning: Read Safety Information

Unique device identifier

Prescription use only

Indicates the authorized  
representative in Switzerland

Formlabs Ohio Inc.  
27800 Lemoyne Rd,  
Millbury, OH 43447  
USA  
+1 617 855 0762

EU Importer  
Formlabs GmbH  
Nalepastr. 18  
12459 Berlin, Germany  
+49 1573 5993322

Curaden AG  
Amlehnstrasse 22  
6010 Kriens  
Switzerland

Before printing, refer to the Manufacturing Guide and dental applications guide for this material at [support.formlabs.com](https://support.formlabs.com) for printing and post-processing recommendations and requirements to ensure the correct and safe usage of this material. For detailed safety and environmental information, the Safety Data Sheet is available at [support.formlabs.com](https://support.formlabs.com).

## I. Introduction And Indications For Use

### INTRODUCTION AND INTENDED PURPOSE

IBT Flex Resin is light-curable polymer-based resin designed for the additive fabrication of biocompatible, short-term use, removable dental appliances such as indirect bonding trays and trays for guided restorative techniques.

### INDICATIONS FOR USE

IBT Flex Resin is a photopolymer resin that can be utilized to 3D print dental appliances such as indirect bonding trays to aid in the placement of brackets/attachments in the treatment of occlusal disharmonies and trays for guided restorative techniques to aid in the molding of composite in the treatment of aesthetic or functional disharmonies.

### CONTRAINDICATIONS

IBT Flex Resin should not be used with known allergy or hypersensitivity history to photopolymer resins. In case of an allergic reaction, stop use immediately, and consult a physician.

### INTENDED PATIENT POPULATION

All patients and populations in which a clinician has prescribed a treatment associated with the intended purpose.

### INTENDED USERS

Dental professionals who fabricate dental appliances.

### CLINICAL BENEFIT

The clinical benefit of IBT Flex Resin is the ability to print patient-specific dental appliances such as indirect bonding trays and trays for guided restorative techniques.

## II. Safety

### RESIN ASSOCIATED WARNINGS

IBT Flex Resin (uncured) contains polymerizable monomers.

- Please read IBT Flex Resin SDS at [support.formlabs.com](https://support.formlabs.com) for hazards associated with liquid resin.

### RESIN USE ASSOCIATED PRECAUTIONS

- Protective glasses and nitrile gloves should be worn while handling resin. For information about the handling of IBT Flex Resin please read Safety Data Sheet at [support.formlabs.com](https://support.formlabs.com).

### INFORMATION FOR PATIENT SAFETY

- The patient should inform the clinician of any known allergy or hypersensitivity history to photopolymer resins.
- In case of an allergic reaction, stop use immediately and consult a medical physician.
- Do not use appliances if damage or cracks are present.
- If signs of extreme wear, deterioration, and/or perforations appear, the appliances should be replaced to avoid material failure and/or potential bodily injury.

## III. Performance Characteristics

For detailed information on the mechanical properties of this material, technical data sheets can be found at [support.formlabs.com](https://support.formlabs.com). IBT Flex Resin has been evaluated in accordance with ISO 10993-1:2020 and ISO 7405:2018, and meets all the biocompatibility requirements for its indicated uses.

## IV. Specific Manufacturing Considerations

Refer to the Manufacturing Guide at [support.formlabs.com](https://support.formlabs.com) for printing recommendations and requirements to ensure the correct and safe usage of this material.

- Hardware: Formlabs 3D Printers and accessories

- Software: Formlabs PreForm
- Part Orientation: Design guidelines for part orientation during printing can be found in the Manufacturing Guide.
- Part Thickness: Design guidelines for part thickness can be found in the Manufacturing Guide.
- Recommended Post-Processing Equipment:
  - Formlabs validated accessories
  - Formlabs validated washing units with approved solvent
  - Formlabs validated curing units

## REQUIREMENTS

For biocompatibility compliance, IBT Flex Resin was validated with a dedicated resin tank, build platform, wash unit, and post-processing equipment that were not mixed with any other resins.

### A. PRINTING

Printing guidelines and validated equipment can be found in the Manufacturing Guide at [support.formlabs.com](https://support.formlabs.com).

### B. WASHING

Washing guidelines and validated equipment can be found in the Manufacturing Guide.

### C. POST-CURING

Post-Curing guidelines and validated equipment can be found in the Manufacturing Guide.

### D. SUPPORT REMOVAL

Support marks can cause abrasion if not removed and polished. See the Manufacturing Guide for polishing recommendations.

### E. CLEANING & DISINFECTION

1. Fully post-processed parts can be cleaned using a dedicated soft toothbrush with neutral soap and room-temperature water. Do not use toothpaste or any other cleaning products on parts printed with IBT Flex Resin. After cleaning, always inspect parts for any cracks. Discard if any damage or cracks are detected.
2. The appliances may be cleaned and disinfected according to facility protocols. Tested disinfection method: soaking the finished appliance in fresh 70% IPA for 5 minutes. Do not leave the part in the alcohol solution for longer than 5 minutes. After disinfection, inspect the part for damage or cracks to ensure that the integrity of the designed part meets performance requirements.

### F. STORAGE

1. Printed appliances should be stored in closed, opaque, or amber containers. Excess light exposure over time may affect the color and performance of printed appliances.
2. Store cartridges at 10°C - 25°C (50°F - 77°F). Do not exceed 25°C (77°F) when in storage.
3. Keep cartridges closed & away from ignition sources.

### G. DISPOSAL

1. Cured resin is non-hazardous and may be disposed of as regular waste.
2. Follow facility protocols for waste that may be considered biohazardous.
3. Unused liquid resin & solvents should be disposed of in accordance with local regulations.
4. Cartridge & contaminated packaging should be disposed of in accordance with local regulations.

### H. ADVERSE EVENTS

To report any serious injury (temporary or serious permanent deterioration of health, a serious public health threat, or death), contact Formlabs Inc. (+1 617 855 0762) or your local Competent Authority.

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