INSTRUCTIONS FOR USE

Surgical Guide Resin

Prepared: 03/22/2021 PRNT-0001 Rev 02

Instructions for use

The following instructions for use are for Formlabs photopolymer Surgical Guide Resin. Basic information about safety and environmental concerns are also included. For more detailed safety and environment information please refer to the Safety Data Sheet, available at dental.formlabs.com. For further information regarding the use of this material, please contact Formlabs.

Prepared: 03/22/2021 PRNT-0001 Rev 02

: Keep away from sunlight

: Consult Instructions for Use

: Batch Code

: Manufacturer

: Authorized representative in the European Community

: Environmental hazard

: Use-by date

: European Conformity

: Catalog Number

: Temperature Limit

: Medical Device

: Warning: Read Safety Information

Formlabs Ohio Inc.
27800 Lemoyne Rd,
Millbury, OH 43477 USA
+1 617 855 0762

EU Importer

Formlabs GmbH
Funkhaus Berlin, Nalepastr. 18
12459 Berlin, Germany
+49 1573 5993322

1. Introduction And Indications For Use

Introduction

Surgical Guide Resin is a light-curable polymer-based resin designed for the fabrication of biocompatible, short term use, removable dental appliances by additive manufacturing.

INTENDED PURPOSE

It is intended for 3D printing endosseous dental implant accessories such as dental surgical guides.

INDICATIONS FOR USE

Surgical Guide Resin is a photopolymer resin indicated for 3D printing biocompatible dental surgical guides to aid in endosseous implant placement surgeries.

CONTRAINDICATIONS

Surgical Guide Resin should not be used in case of any known allergy or hypersensitivity history to photopolymer resin materials. In case of an allergic reaction, stop use immediately, and consult a medical physician.

INTENDED PATIENT POPULATION

All patients and populations in which a clinician has prescribed a treatment associated with the intended purpose.

INTENDED USERS

Dental professionals who fabricate dental appliances.

CLINICAL BENEFIT

Clinical benefits may vary depending on patient case and appliance design:

Appliance	Clinical Benefit
Surgical Guide	Accuracy: Assists dentists or oral surgeons to place or direct surgical dental implants into the bone with increased accuracy and precision. Operative Time: Reduction of chair time, and invasiveness of dental implant surgeries.

2. Safety

RESIN ASSOCIATED WARNINGS

Surgical Guide Resin (uncured) contains polymerizable monomers which may cause skin irritation or other allergic reactions in susceptible persons. If resin contacts skin, wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.

1. **Eye contact:** High vapor concentration may cause irritation.
2. **Skin contact:** May cause sensitization by skin contact. Irritating to skin. Repeated and/or prolonged contact may cause dermatitis.
3. **Inhalation:** Irritating to the respiratory system. Prolonged or repeated exposure may cause: headache, drowsiness, nausea, weakness (severity of effects depends on the extent of exposure).
4. **Ingestion:** Low oral toxicity, but ingestion may cause irritation of the gastrointestinal tract.
5. **Protection:** Protective glasses and nitrile gloves should be worn while handling Surgical Guide Resin. Detailed information about the handling of Surgical Guide Resin can be found in the Safety Data Sheets at dental.formlabs.com.

RESIN ASSOCIATED PRECAUTIONS

1. When washing the printed part with solvent, it should be in a properly ventilated environment with proper protective masks and gloves.
2. Expired or unused Surgical Guide Resin shall be disposed in accordance with local regulations.
3. IPA shall be disposed of in accordance with local regulations
4. Disposal recommendations in Section 5. G.

INFORMATION FOR PATIENT SAFETY

1. The patient should inform the clinician of any known allergy or hypersensitivity history to photopolymer resin materials.
2. In case of an allergic reaction, stop use immediately, and consult a medical physician.

3. Hardware And Printing Parameters

- a. **Hardware:** Formlabs SLA 3D Printer (Laser Wavelength: 405 nm)
- b. **Software:** Formlabs PreForm
 - STL file import
 - Manual / Automatic rotation and placement
 - Manual / Automatic generation of supports
- c. **Printing Parameters**
 - Layer thickness: 50um and 100 µm
 - Optimal orientation: Horizontal orientation with the intaglio surface facing away from the build platform
 - Minimum thickness of 2 mm
- d. **Recommended Post-Processing Equipment:**
 - Formlabs Form Wash
 - Isopropyl alcohol (IPA) ≥ 99%
 - Formlabs Form Cure

4. Performance Characteristics

Mechanical Properties	Post Cured Results ¹	Method
Flexural Strength	> 102 MPa	ASTM D638
Flexural Modulus	> 2400 MPa	ASTM D790

Surgical Guide Resin has been evaluated in accordance with EN ISO 10993-1:2020, *Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process*, and EN ISO 7405:2018, *Dentistry - Evaluation of biocompatibility of medical devices used in dentistry*, and passed the requirements for the following biocompatibility risks:

¹ Material properties may vary based on part geometry, print orientation, print settings, temperature, and disinfection or sterilization methods used.

ISO Standard	Description ²
EN ISO 10993-5	Not cytotoxic
EN ISO 10993-10	Not an irritant
EN ISO 10993-10	Not a sensitizer

5. Specific Manufacturing Considerations

NOTIFICATION

The device specifications have been validated using the printer process parameters indicated below.

REQUIREMENTS

Use dedicated accessories for Surgical Guide Resin. For biocompatibility compliance, Surgical Guide Resin requires a dedicated resin tank, build platform, Form Wash, and finishing kit, which should not be mixed with any other resins.

A. PRINTING

1. **Shake cartridge:** Prior to use, shake cartridge well. Color deviations and print failures may occur if the cartridge is shaken insufficiently.
2. **Set up:** Insert resin cartridge into a compatible Formlabs 3D printer.
3. **Printing:**
 - a. Prepare a print job using PreForm software. Import desired part STL file. Orient and generate supports. For recommendations on print orientation and support placement, visit dental.formlabs.com.
 - b. Send print job to printer. Begin print by selecting a print job from the print menu. Follow any prompts or dialogs shown on the printer screen. Printer will automatically complete the print.
4. **Cartridge Reuse:**
 - a. After the print is completed, store unused resin cartridges as instructed in Section F below.
 - b. For a new print job, make sure to shake the cartridge well and follow all instructions as recommended.

² Surgical Guide Resin was tested at NAMSA World Headquarters, OH, USA.

B. POST-PROCESSING

Warning: During production, inspect the printed appliance for cracks, rough edges, or any other artifacts that may affect finished device performance or patient comfort and safety. Part inspection is recommended after every fabrication step, for example but not limited to, after printing, support removal, cleaning, and/or sterilization when applicable.

1. **Part removal:**

- a. Remove the build platform from the printer.
- b. Printed parts can be removed from the build platform before or after cleaning in a Form Wash. To remove, wedge the part removal tool under the printed part raft, and rotate the tool. For detailed techniques, visit support.formlabs.com.

2. **Rinsing:** Required - Place the printed parts in a Form Wash filled with Isopropyl Alcohol (IPA, $\geq 99\%$) and wash for 20 minutes or until clean.

Note: Printed parts with complex geometries may require additional washing.

3. **Drying:**

- a. Remove parts from IPA and leave to air dry at room temperature for **at least** 30 minutes. **Note:** Dry times can vary depending on design of parts and ambient conditions.
- b. Inspect printed parts to ensure that parts are clean and dry. No residual alcohol, excess liquid resin or residue particles must remain on surface before proceeding to subsequent steps.
- c. If residual alcohol is still present, dry parts longer. If resin residue is still visible, rewash parts until clean and dry.

4. **Post cure:** Required - Place the dried printed parts in a Form Cure and post cure as recommended below:

Printer	Temperature	Time
Form 2	60°C	30 minutes
Form 3B	70°C	30 minutes
Form 3BL	70°C	30 minutes

Allow Form Cure to cool down to room temperature between cure cycles. **Note:** During post-curing, a color change from translucent yellow to translucent orange will occur.

5. **Support removal:**

- a. Remove supports using a cutting disk and handpiece, or with other part removal tools.
- b. Inspect the printed parts for any cracks. Discard if any damage or cracks are detected.

C. ASSEMBLY

1. Press the metal surgical drill guide sleeve(s) in the corresponding hole(s) of the printed surgical guide.
2. Ensure to insert metal sleeves that are validated and deemed compatible for use.

D. POLISHING

1. If any rough marks are left on the printed surgical guide surface after support removal, smooth down the support surfaces with a bur and handpiece to improve patient comfort.
2. Polish the printed surgical guides using typical dental polishing methods prior to patient use.

E. DISINFECTION AND STERILIZATION

1. The surgical guide may be cleaned, disinfected and sterilized according to facility protocols. Tested methods of disinfection include: soaking the finished surgical guide in fresh 70% IPA for 5 minutes. **Note:** Do not leave the part in an alcohol solution for an extended period as this may adversely affect performance.
2. Surgical guides may be steam sterilized according to CDC recommended cycles (132°C / 270°F for 4 minutes in a pre-vacuum steam sterilizer or 30 minutes at 121°C / 250°F in a gravity displacement autoclave), or according to facility/autoclave manufacturer's protocol, as long as cycles do not exceed 20 minutes for 134°C / 273°F or 30 minutes at 121°C / 250°F.

3. Autoclave cycles should include a dry cycle to best maintain accuracy. For example, wrapped instruments sterilized in a prevacuum autoclave should be dried for 20-30 minutes according to CDC recommendations.³ Longer or hotter autoclave cycles than those listed above may result in degradation of physical properties and accuracy. **Note:** A color shift will be observed after autoclaving, this is normal.
4. After disinfection and sterilization, inspect the printed parts for cracks to ensure the integrity of the surgical guide.

F. STORAGE

1. When not in use, place printed parts in closed, opaque or amber containers.
2. Store in a cool, dry place out of direct sunlight. Excess light exposure overtime may affect the performance of the printed parts.
3. Store the cartridges at 10°C - 25°C (50°F - 77°F).
4. Do not exceed 25°C (77°F) when in storage.
5. Keep away from ignition sources.

G. DISPOSAL

1. Any cured resin is non-hazardous and may be disposed of as regular waste.
2. Follow clinic or facility protocols for waste that may be considered biohazardous.
3. Liquid resin should be disposed of in accordance with government regulations (community, regional, national).
 - a. Contact a licensed professional waste disposal service to dispose of liquid resin.
 - b. Do not allow waste to enter storm or sewer drainage systems.
 - c. Avoid release into the environment.
 - d. Contaminated packaging: Dispose of cartridge or any unused product per local waste management guidelines.

³ Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008, Centers for Disease Control and Prevention (CDC), Updated: May 2019

H. ADVERSE EVENTS

To report any serious injury (temporary or permanent serious deterioration of health, a serious public health threat, or death), contact Formlabs Inc. (+1 617 855 0762) or your local Competent Authority ⁴.

⁴ Local EU Competent Authorities listed here:

https://ec.europa.eu/health/sites/health/files/md_sector/docs/md_vigilance_contact_points.pdf

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